

Datasheet: PPP019

Description:	RECOMBINANT PIG CXCL9
Name:	CXCL9
Other names:	MIG
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			
Western Blotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified recombinant protein - lyophilised
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Preservative Stabilisers	10% Trehalose
External Database Links	UniProt: B0FYK2 Related reagents

Product Information	Recombinant Pig CXCL9 is a recombinant protein corresponding to aa 23 - 126 of porcine CXCL9, also known as MIG (monokine induced by interferon-gamma). CXCL9 is a member of the CXC family of chemokines, secreted by interferon γ stimulated monocytes, macrophages and endothelial cells, that signals through binding to CD183 (CXCR3). CXCL9 acts as a Th1 (type 1 helper T) cell chemoattractant and plays a role in the growth, activation and movement of cells associated with immune and inflammatory responses, and in tumor growth inhibition.
----------------------------	---

Protein Molecular	12.1 kDa
--------------------------	----------

Weight

Purity >95% by SDS PAGE analysis

ELISA This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2387](#) as a coating antibody and [AHP2387B](#) as a detection antibody.

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted.

Shelf Life 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10591 available at:
10591: <https://www.bio-rad-antibodies.com/uploads/MSDS/10591.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M287056:160523'

Printed on 02 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)