

Datasheet: PPP018

Description:	RECOMBINANT PIG MIP-1 BETA
Name:	MIP-1 BETA
Other names:	CCL4
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			
Western Blotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose
External Database Links	UniProt: Q711P4 Related reagents Entrez Gene: 396668 CCL4 Related reagents
Synonyms	MIP1B, SCY4A
Product Information	Recombinant Pig MIP-1 beta (macrophage inflammatory protein-1 beta), also known as C-C motif

chemokine 4 (CCL4), is a recombinant protein corresponding to aa 24-92 of native pig MIP-1 β . Porcine MIP-1 β is a secreted protein belonging to the intercrine beta (chemokine CC) family.

Porcine MIP-1 β functions in inflammatory responses and as a chemoattractant for various cell types including NK cells and monocytes.

Protein Molecular Weight	9.3 kDa
Purity	>95% by SDS PAGE analysis
ELISA	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with AHP2386 as a coating antibody and AHP2386B as a detection antibody.
Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10591 available at: 10591: https://www.bio-rad-antibodies.com/uploads/MSDS/10591.pdf
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M287053:160523'

Printed on 02 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)