

Datasheet: PRABP008

Description:	RECOMBINANT RABBIT MCP-1
Name:	MCP-1
Other names:	CCL2
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			
Western Blotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose
External Database Links	UniProt: P28292 Related reagents Entrez Gene: 100009130 CCL2 Related reagents
Synonyms	SCYA2
Product Information	Recombinant Rabbit MCP-1 is an 11.4 kDa peptide corresponding to aa 24 – 125 of native rabbit

monocyte chemotactic protein 1 (MCP-1). Also known as CCL2 , MCP-1 is a secreted chemokine involved in the recruitment of monocytes to sites of injury or inflammation. MCP-1 plays an important role in the inflammatory response of blood monocytes and tissue macrophages.

Protein Molecular Weight	11.4 kDa
Purity	>95% by SDS PAGE analysis
ELISA	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with AHP2994 as a coating antibody and AHP2994B as a detection antibody.
Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10591 available at: 10591: https://www.bio-rad-antibodies.com/uploads/MSDS/10591.pdf
Regulatory	For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M294997:161012'

Printed on 02 May 2018