

Datasheet: STAR153

Description:	GOAT F(ab') ₂ IgG NEGATIVE CONTROL
Specificity:	GOAT F(ab') ₂ IgG NEGATIVE CONTROL
Format:	Purified
Product Type:	Negative/Isotype Control
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Negative Control
Product Form	Purified IgG - liquid
Buffer Solution	Borate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	F(ab') ₂ concentration 5mg/ml
Specificity	Goat F(ab')₂ negative control reagent is intended as a negative control for use in flow cytometry and ELISA to assess potential non-specific binding by goat F(ab') ₂ IgG antibodies.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch.
Health And Safety	Material Safety Datasheet documentation available at:

Information

Material Safety Datasheet Documentation #10077 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10077.pdf>

Regulatory

For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M310317:170915'

Printed on 15 Sep 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)