

Datasheet: STAR133A

Description:	GOAT ANTI MOUSE IgG2a:Alk. Phos.
Specificity:	IgG2a
Format:	Alk. Phos.
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Immunohistology - Frozen	■			
Immunohistology - Paraffin	■			
ELISA	■			1/2000 - 1/4000
Western Blotting	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Species Cross Reactivity	Does not react with:Human
Product Form	Purified IgG conjugated to Alkaline Phosphatase - liquid
Antiserum Preparation	Antisera to mouse IgG2a were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	TRIS buffered saline, 1mM MgCl ₂
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃) 50% Glycerol
Immunogen	IgG2a paraproteins from BALB/c mice

External Database Links

UniProt:

P01865	Related reagents
P01863	Related reagents
P01864	Related reagents

Entrez Gene:

[380793](#) Igh-1a [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)

Specificity **Goat anti Mouse IgG2a antibody** recognizes Mouse IgG2a. This antibody has been cross absorbed against mouse IgM, IgG1, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG2a antibody shows minimal cross-reactivity with human immunoglobulins.

References

1. Knipping, K. *et al.* (2011) A gastrointestinal rotavirus infection mouse model for immune modulation studies. [Virol J. 8: 109.](#)
2. Bagai, U. & Pawar, A. (2013) A blood stage fraction of Plasmodium berghei induces protective and long lasting immune response in BALB/c mice. [Parasitol Int. 62 \(3\): 329-36.](#)
3. Hwang, S.R. *et al.* (2015) Altered expression levels of neurodevelopmental proteins in fetal brains of BTBR T+tf/J mice with autism-like behavioral characteristics. [J Toxicol Environ Health A. 78 \(8\): 516-23.](#)
4. Zhao, Z. *et al.* (2015) Multiple B-cell epitope vaccine induces a *Staphylococcus* enterotoxin B-specific IgG1 protective response against MRSA infection. [Sci Rep. 5: 12371.](#)

Storage Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10322 available at:
10322: <https://www.bio-rad-antibodies.com/uploads/MSDS/10322.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M314552:180412'

Printed on 02 May 2018