

Datasheet: STAR133PE

Description:	GOAT ANTI MOUSE IgG2a:RPE
Specificity:	IgG2a
Format:	RPE
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/500

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse		
Species Cross Reactivity	Does not react with:Human		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578

Antiserum Preparation Antisera to mouse IgG2a were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.

Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide Sucrose
Approx. Protein Concentrations	IgG concentration 0.5mg/ml
Immunogen	IgG2a paraproteins from BALB/c mice
External Database Links	UniProt:

[P01865](#) [Related reagents](#)

[P01863](#) [Related reagents](#)

[P01864](#) [Related reagents](#)

Entrez Gene:

[380793](#) Igh-1a [Related reagents](#)

[380793](#) Igh-1a [Related reagents](#)

[380793](#) Igh-1a [Related reagents](#)

Specificity	Goat anti Mouse IgG2a antibody recognizes Mouse IgG2a. This antibody has been cross absorbed against mouse IgM, IgG1, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG2a antibody shows minimal cross-reactivity with human immunoglobulins.
--------------------	--

Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
-----------------------	---

References	1. Knipping, K. <i>et al.</i> (2011) A gastrointestinal rotavirus infection mouse model for immune modulation studies. Virol J. 8: 109. 2. Bagai, U. & Pawar, A. (2013) A blood stage fraction of Plasmodium berghei induces protective and long lasting immune response in BALB/c mice. Parasitol Int. 62 (3): 329-36. 3. Hwang, S.R. <i>et al.</i> (2015) Altered expression levels of neurodevelopmental proteins in fetal brains of BTBR T+tf/J mice with autism-like behavioral characteristics. J Toxicol Environ Health A. 78 (8): 516-23. 4. Zhao, Z. <i>et al.</i> (2015) Multiple B-cell epitope vaccine induces a <i>Staphylococcus</i> enterotoxin B-specific IgG1 protective response against MRSA infection. Sci Rep. 5: 12371.
-------------------	---

Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

Shelf Life	6 months from date of despatch.
-------------------	---------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10045 available at: 10045: https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M314578:180412'

Printed on 02 May 2018