

## Datasheet: STAR116

|                      |                     |
|----------------------|---------------------|
| <b>Description:</b>  | GOAT ANTI RAT IgM   |
| <b>Specificity:</b>  | IgM                 |
| <b>Format:</b>       | Purified            |
| <b>Product Type:</b> | Polyclonal Antibody |
| <b>Isotype:</b>      | Polyclonal IgG      |
| <b>Quantity:</b>     | 0.5 mg              |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                            | Yes | No | Not Determined | Suggested Dilution |
|----------------------------|-----|----|----------------|--------------------|
| Flow Cytometry             | ■   |    |                |                    |
| Immunohistology - Frozen   |     |    | ■              |                    |
| Immunohistology - Paraffin |     |    | ■              |                    |
| ELISA                      | ■   |    |                | 1ug/ml - 10ug/ml   |
| Immunoprecipitation        |     |    | ■              |                    |
| Western Blotting           | ■   |    |                |                    |

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

|                                       |                                                                                                                                                                                   |
|---------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Target Species</b>                 | Rat                                                                                                                                                                               |
| <b>Product Form</b>                   | Purified IgG - liquid                                                                                                                                                             |
| <b>Antiserum Preparation</b>          | Antisera to rat IgM were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.           |
| <b>Buffer Solution</b>                | Phosphate buffered saline                                                                                                                                                         |
| <b>Preservative Stabilisers</b>       | 0.09% Sodium Azide                                                                                                                                                                |
| <b>Approx. Protein Concentrations</b> | IgG concentration 1.0 mg/ml                                                                                                                                                       |
| <b>Specificity</b>                    | <b>Goat anti Rat IgM antibody</b> recognizes rat IgMs assessed by immunoelectrophoresis and ELISA.<br><br>Goat anti Rat IgM antibody may cross react with IgM from other species. |
| <b>Flow Cytometry</b>                 | Use 10ul of the suggested working dilution to stain $1 \times 10^6$ cells in 100ul.                                                                                               |

## References

1. Terrell, S.P. *et al.* (2012) Glomerulonephropathy in aged captive Key Largo woodrats (*Neotoma floridana smalli*). [Vet Pathol. 49 \(4\): 710-6.](#)
2. Kesik-brodacka, M. *et al.* (2017) Immune response of rats vaccinated orally with various plant-expressed recombinant cysteine proteinase constructs when challenged with *Fasciola hepatica* metacercariae. [PLoS Negl Trop Dis. 11 \(3\): e0005451.](#)
3. Bouma, H.R. *et al.* (2013) Hibernation is associated with depression of T-cell independent humoral immune responses in the 13-lined ground squirrel. [Dev Comp Immunol. 39 \(3\): 154-60.](#)

---

## Storage

Store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

---

## Shelf Life

12 months from date of despatch.

---

## Health And Safety Information

Material Safety Datasheet documentation #10040 available at:  
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

---

## Regulatory

For research purposes only

---

## Related Products

### Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

#### North & South America

Tel: +1 800 265 7376  
Fax: +1 919 878 3751  
Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

#### Worldwide

Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739  
Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

#### Europe

Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50  
Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

'M316896:180612'

Printed on 20 Jun 2018