

Datasheet: STAR79B

Description:	GOAT ANTI HAMSTER IgG:Biotin
Specificity:	IgG
Format:	Biotin
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.8 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/25 - 1/50
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			1/1000 - 1/5000

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Hamster
Product Form	Purified IgG conjugated to Biotin - liquid
Antiserum Preparation	Antisera to hamster IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Stabilisers	1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.8 mg/ml
Immunogen	Hamster IgG.
Specificity	Goat anti Hamster IgG antibody recognizes hamster IgG (H+L) and has been adsorbed against both mouse and rat immunoglobulins to minimise cross-reactivity.
Flow Cytometry	Use 50ul of the suggested working dilution to label 10^6 cells in 100ul.

References

1. Samant, M. *et al.* (2009) Immunization with the DNA-encoding N-terminal domain of proteophosphoglycan of *Leishmania donovani* generates Th1-type immunoprotective response against experimental visceral leishmaniasis. [J Immunol. 183: 470-9.](#)
2. Moore, G.T. *et al.* (2008) Glycosylation changes in hFUT1 transgenic mice increase TCR signaling and apoptosis resulting in thymocyte maturation arrest. [Mol Immunol. 45: 2401-10.](#)
3. Forster, K.M. *et al.* (2015) DNA prime-protein boost based vaccination with a conserved region of leptospiral immunoglobulin-like A and B proteins enhances protection against leptospirosis. [Mem Inst Oswaldo Cruz. 110 \(8\): 989-95.](#)
4. Verma R *et al.* (2015) Cross reactive molecules of human lymphatic filaria *Brugia malayi* inhibit *Leishmania donovani* infection in hamsters. [Acta Trop. 152: 103-11.](#)
5. Wiśniewski, M. *et al.* (2016) Hamsters vaccinated with Ace-mep-7 DNA vaccine produced protective immunity against *Ancylostoma ceylanicum* infection. [Exp Parasitol. 163: 1-7.](#)
6. Bacelo, K.L. *et al.* (2014) Xanthan gum as an adjuvant in a subunit vaccine preparation against leptospirosis. [Biomed Res Int. 2014: 636491.](#)

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

Health And Safety Information

Material Safety Datasheet documentation #10041 available at:
10041: <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory

For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M316888:180612'

Printed on 21 Jun 2018