

## Datasheet: STAR97PE

<b>Description:</b>	GOAT F(ab') <sub>2</sub> ANTI HUMAN IgG:RPE
<b>Specificity:</b>	IgG (H/L)
<b>Format:</b>	RPE
<b>Product Type:</b>	Polyclonal Antibody
<b>Isotype:</b>	Polyclonal IgG
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human		
<b>Product Form</b>	F(ab') <sub>2</sub> fragment of IgG conjugated to R. Phycoerythrin (RPE) - lyophilised		
<b>Reconstitution</b>	Reconstitute with 1 ml distilled water		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	RPE 488nm laser	496	578

**Antiserum Preparation** Antisera to human IgG were raised by repeated immunisation of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography. F(ab')<sub>2</sub> fragments were prepared by pepsin digestion.

<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide
<b>Stabilisers</b>	1% Bovine Serum Albumin 5% Sucrose
<b>Immunogen</b>	Human IgG

### External Database Links

**UniProt:**  
[P01857](#)    [Related reagents](#)  
[P01834](#)    [Related reagents](#)

[P01860](#)    [Related reagents](#)  
[P01859](#)    [Related reagents](#)  
[P01861](#)    [Related reagents](#)  
[P0CG04](#)    [Related reagents](#)

**Entrez Gene:**

[3500](#)    IGHG1    [Related reagents](#)  
[3514](#)    IGKC    [Related reagents](#)  
[3501](#)    IGHG2    [Related reagents](#)  
[3502](#)    IGHG3    [Related reagents](#)  
[3503](#)    IGHG4    [Related reagents](#)  
[28815](#)    IGLV2-14    [Related reagents](#)

---

**Specificity**

**Goat F(ab')<sub>2</sub> anti Human IgG antibody** recognizes human IgG and demonstrates <1% cross-reactivity with other human immunoglobulin classes or subclasses.

Goat F(ab')<sub>2</sub> anti Human IgG antibody does not cross-react with Goat IgG F(c) or mouse, rat, bovine and horse serum proteins, as determined by immunoelectrophoresis.

Goat F(ab')<sub>2</sub> anti Human IgG antibody has been used for the detection of the antibody response to an endogenous retrovirus GAG matrix in patients with rheumatoid arthritis by ELISA ([Nelson \*et al.\* 2014](#)).

---

**Flow Cytometry**

Use 50ul of the suggested working dilution to label 10<sup>6</sup> cells in 100ul.

---

**References**

1. Nelson, P.N. *et al.* (2014) Rheumatoid Arthritis is Associated with IgG Antibodies to Human Endogenous Retrovirus Gag Matrix: A Potential Pathogenic Mechanism of Disease? [J Rheumatol. 41: 1952-60.](#)

---

**Storage**

Store at +4°C.

DO NOT FREEZE

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Shelf Life**

12 months from date of reconstitution.

---

**Health And Safety Information**

Material Safety Datasheet documentation #10075 available at:  
10075: <https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf>

---

**Regulatory**

For research purposes only

---

## Related Products

### Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

**North & South America**    Tel: +1 800 265 7376  
Fax: +1 919 878 3751

**Worldwide**    Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

**Europe**    Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

'M314452:180412'

**Printed on 20 Jun 2018**

---