

Datasheet: VMA00138

Description:	MOUSE ANTI CD80
Specificity:	CD80
Format:	Purified
Product Type:	PrecisionAb™ Monoclonal
Clone:	OT12B11
Isotype:	IgG2a
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Mouse monoclonal antibody purified by affinity chromatography from ascites.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 50% Glycerol
Immunogen	Full length recombinant human CD80 (NP_005182) produced in HEK293T cells
External Database Links	UniProt: P33681 Related reagents Entrez Gene: 941 CD80 Related reagents
Synonyms	CD28LG, CD28LG1, LAB7

Specificity	Mouse anti Human CD80 antibody recognizes CD80, also known as B-lymphocyte activation antigen B7, CD80 antigen (CD28 antigen ligand 1, B7-1 antigen), CTLA-4 counter-receptor B7.1, T-lymphocyte activation antigen CD80, activation B7-1 antigen, costimulatory factor CD80 and costimulatory molecule variant IgV-CD80. CD80 is a membrane receptor that is activated by the binding of CD28 or CTLA-4. The activated protein induces T-cell proliferation and cytokine production. This protein can act as a receptor for adenovirus subgroup B and may play a role in lupus neuropathy (provided by RefSeq, Aug 2011). Mouse anti Human CD80 antibody detects a band of 60 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.
Western Blotting	Anti CD80 detects a band of approximately 60 kDa in Raji cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10048 available at: Antibody (10048): https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M284501:160321'

Printed on 30 Apr 2018