

Datasheet: VMA00009

Description:	HUMAN ANTI CYTOKERATIN 19/20
Specificity:	CYTOKERATIN 19/20
Format:	Purified
Product Type:	PrecisionAb™ Monoclonal
Isotype:	HuCAL Fab bivalent
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	A bivalent human recombinant Fab selected from the HuCAL® GOLD phage display library, expressed in <i>E. coli</i> . This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This product is supplied as a liquid.
Preparation	Recombinant human monoclonal antibody purified by metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Recombinant human cytokерatin 20
External Database Links	<p>UniProt:</p> <p>P35900 Related reagents</p> <p>P08727 Related reagents</p> <p>Entrez Gene:</p> <p>54474 KRT20 Related reagents</p> <p>3880 KRT19 Related reagents</p>

Specificity	Human anti Human cytokeratin 19/20 antibody recognizes human cytokeratin 20, and has been demonstrated to cross-react with cytokeratin 19. In immunohistology reactivity is seen both with tissues expressing only CK19 and with those expressing only CK20. Human anti Human cytokeratin 19/20 antibody recognizes cytokeratins 19 and 20 as bands of ~40 kDa and ~45 kDa respectively by western blotting under reducing conditions.
Western Blotting	Anti cytokeratin 19/20 antibody detects bands of approximately 40 kDa and 45 kDa in HeLa cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. PrecisionAb™ is a trademark of Bio-Rad Laboratories. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Goat Anti Human IgG F(ab')₂ (STAR126...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M311371:171113'

Printed on 02 May 2018