

## Datasheet: VMA00147

<b>Description:</b>	MOUSE ANTI ErbB2
<b>Specificity:</b>	ErbB2
<b>Format:</b>	Purified
<b>Product Type:</b>	PrecisionAb™ Monoclonal
<b>Clone:</b>	OT14F10
<b>Isotype:</b>	IgG2a
<b>Quantity:</b>	100 µl

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

**PrecisionAb antibodies have been extensively [validated for the western blot application](#).** The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

<b>Target Species</b>	Human
<b>Species Cross Reactivity</b>	Reacts with: Mouse, Rat <b>N.B.</b> Antibody reactivity and working conditions may vary between species.
<b>Product Form</b>	Purified IgG - liquid
<b>Preparation</b>	Mouse monoclonal antibody purified by affinity chromatography from ascites.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.09% Sodium Azide (NaN <sub>3</sub> ) 1% Bovine Serum Albumin 50% Glycerol
<b>Immunogen</b>	Recombinant protein corresponding to amino acids 676-1255 of human ErbB2 (NP_004439) produced in HEK293T cells
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P04626</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b>

---

<b>Synonyms</b>	HER2, MLN19, NEU, NGL
-----------------	-----------------------

---

<b>Specificity</b>	<p><b>Mouse anti Human ErbB2 antibody</b> recognizes ErbB2, also known as c-erb B2/neu protein, herstatin, metastatic lymph node gene 19 protein, neuro/glioblastoma derived oncogene homolog, neuroblastoma/glioblastoma derived oncogene homolog, p185erbB2, proto-oncogene Neu, proto-oncogene c-ErbB-2, receptor tyrosine-protein kinase erbB-2, tyrosine kinase-type cell surface receptor HER2, v-erb-b2 avian erythroblastic leukemia viral oncoprotein 2, v-erb-b2 erythroblastic leukemia viral oncogene homolog 2 and neuro/glioblastoma derived oncogene homolog.</p> <p>Encoded by the ERBB2 gene, ErbB2 is a member of the epidermal growth factor (EGF) receptor family of receptor tyrosine kinases. This protein has no ligand binding domain of its own and therefore cannot bind growth factors. However, it does bind tightly to other ligand-bound EGF receptor family members to form a heterodimer, stabilizing ligand binding and enhancing kinase-mediated activation of downstream signalling pathways, such as those involving mitogen-activated protein kinase and phosphatidylinositol-3 kinase. Allelic variations at amino acid positions 654 and 655 of isoform a (positions 624 and 625 of isoform b) have been reported, with the most common allele, Ile654/Ile655, shown here. Amplification and/or overexpression of ERBB2 has been reported in numerous cancers, including breast and ovarian tumors. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized (provided by RefSeq, Jul 2008).</p> <p>Mouse anti Human ErbB2 antibody detects a band of 185 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
--------------------	--

---

<b>Western Blotting</b>	Anti ErbB2 detects a band of approximately 185 kDa in HepG2 cell lysates.
-------------------------	---

---

<b>Instructions For Use</b>	Please refer to the <a href="#">PrecisionAb western blotting protocol</a> . For additional information on secondary antibody dilution and exposure time see product web page.
-----------------------------	---

---

<b>Storage</b>	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
----------------	---

---

<b>Shelf Life</b>	As supplied, 12 months from date of despatch.
-------------------	---

---

<b>Acknowledgements</b>	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
-------------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10048 available at: Antibody (10048): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10048.pdf</a>
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

### Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA929\)](#)

'M284507:160321'

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

**Printed on 02 May 2018**

---

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)