

Datasheet: VMA00442

Description:	HUMAN ANTI HEAT SHOCK 70 kDa PROTEIN 4
Specificity:	HEAT SHOCK 70 kDa PROTEIN 4
Format:	Purified
Product Type:	PrecisionAb™ Monoclonal
Isotype:	HuCAL Fab bivalent
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	A bivalent human recombinant Fab selected from the HuCAL® phage display library, expressed in <i>E. coli</i> . This Fab fragment is bivalent by dimerization of the bacterial alkaline phosphatase fusion protein. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.
Preparation	Human recombinant monoclonal antibody purified by StrepTactin affinity chromatography from <i>E. coli</i> extracts
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	GST-tagged heat shock 70 kDa protein 4 fusion protein
External Database Links	UniProt: P34932 Related reagents Entrez Gene:

Synonyms	APG2
Specificity	Human anti Human heat shock 70 kDa protein 4 antibody recognizes heat shock 70 kDa protein 4, also known as HSPA4, epididymis secretory sperm binding protein Li 5a and heat shock 70-related protein APG-2. Human anti Human heat shock 70 kDa protein 4 antibody detects a band of 110 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.
Western Blotting	Anti heat shock 70 kDa protein 4 detects a band of approximately 110 kDa in MCF7 cell lysates
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. PrecisionAb™ is a trademark of Bio-Rad Laboratories. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Human IgG F(ab')₂ (STAR126...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M311381:171113'

Printed on 01 May 2018