

Datasheet: VMA00010

Description:	HUMAN ANTI MEK 1
Specificity:	MEK 1
Other names:	Mitogen Activated Protein Kinase 1
Format:	Purified
Product Type:	PrecisionAb™ Monoclonal
Isotype:	HuCAL Fab bivalent
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	A bivalent human recombinant Fab selected from the HuCAL® GOLD phage display library, expressed in <i>E. coli</i> . This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a myc-tag (EQKLISEEDL) and a his-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This product is supplied as a liquid.
Preparation	Recombinant human monoclonal antibody purified by metal chelate affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Full length recombinant human MEK 1 expressed in Sf9 cells
External Database Links	<p>UniProt: Q02750 Related reagents</p> <p>Entrez Gene: 5604 MAP2K1 Related reagents</p>

Synonyms	MEK1, PRKMK1
Specificity	<p>Human anti Human MEK 1 antibody recognizes human dual specificity mitogen-activated protein kinase kinase 1, also known as MEK 1 or MAPKK1. MEK 1 is a 393 amino acid, dual specificity protein kinase which has the ability to phosphorylate both threonine and tyrosine residues. MEK1 acts downstream of Ras and Raf and is responsible for the phosphorylation of extracellular signal-regulated kinase (ERK) mitogen-activated protein kinase (MAPK). It is an integral component of the MAP kinase signalling cascade that regulates cell growth and differentiation, and may also play a role in brain synaptic plasticity.</p> <p>Human anti Human MEK 1 antibody recognizes human MEK1 as a band of ~42 kDa in multiple human cell line lysates by western blotting under reducing conditions.</p>
Western Blotting	Anti MEK 1 antibody detects a band of approximately 42 kDa in Jurkat cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. PrecisionAb™ is a trademark of Bio-Rad Laboratories. His-tag is a registered trademark of EMD Biosciences.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Goat Anti Human IgG F(ab')₂ (STAR126...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M311373:171113'

Printed on 02 May 2018