

Datasheet: VMA00032K

Description:	TRF1 ANTIBODY WITH CONTROL LYSATE
Specificity:	TRF1
Other names:	TERF1
Format:	Purified
Product Type:	PrecisionAb™ Monoclonal
Isotype:	IgG1
Quantity:	10 Westerns

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#).

The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	100µl Mouse monoclonal antibody purified by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Recombinant TRF1
External Database Links	<p>UniProt: P54274 Related reagents</p> <p>Entrez Gene: 7013 TRF1 Related reagents</p>
Synonyms	PIN2, TRBF1, TRF, TRF1

Specificity	<p>Mouse anti Human TRF1 antibody recognizes human telomeric repeat-binding factor 1, also known as TTAGGG repeat binding factor 1 (TRF1) or NIMA-interacting protein 2. TRF1 is a ubiquitously expressed 439 amino acid DNA binding protein and negative regulator of telomere length containing a single Myb domain. Along with tankyrase-1, TIN2, and PINX1, TRF1 binds the double-stranded TTAGGG repeat of telomeres and regulates telomere length. Protection of telomeres-1 (POT1) is a single-stranded telomeric DNA-binding protein which controls the elongation of telomeres by telomerase. TRF1 mediates the loading of POT1 onto single stranded DNA, conveying information about telomere length to the telomere terminus, where telomerase is regulated. If the telomere is long enough, telomerase is inhibited.</p> <p>Additionally, in conjunction with POT1 and TTAGGG repeat binding factor 2 (TRF2), TRF1 recruits factors to the telomere, inducing the formation of the protective t-loop. TRF1 along with TRF2, POT1, TIN2, TPP1 and Rap1 form the shelterin complex which distinguishes telomeres from sites of DNA damage (Wellinger 2009). TRF1 is thought to play a role in leukaemia (Shi et al. 2006).</p> <p>Mouse anti Human TRF1 antibody recognizes human TRF1 as a single band of ~60 kDa in multiple human cell line lysates by western blotting under reducing conditions.</p>
Western Blotting	Anti TRF1 antibody detects a band of approximately 60 kDa in K562 cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Lysate Composition	400µg K562 lysate lyophilized in RIPA buffer.
Lysate Reconstitution	<p>- If using DDT reconstitute the lyophilized lysate with 190µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 10µl 2M DTT.</p> <p>- If using BME reconstitute the lyophilized lysate with 180µl DI H₂O, add 200µl 2x Laemmli Sample Buffer and 20µl BME.</p> <p>Heat at 95°C for 5 minutes. For 10 well mini gels load 25µl. For other gel and comb formats please refer to the PrecisionAb western blotting protocol.</p>
Storage	<p>Antibody: Store undiluted at -20°C, avoiding repeated freeze thaw cycles.</p> <p>Lysate: Store lyophilized lysate at -20°C. After reconstitution aliquot and store at -20°C for up to 3 months or at -80°C for longer term storage.</p>
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	<p>Material Safety Datasheet Documentation #10040 for antibody available at https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf</p> <p>Lysate Material Safety Datasheet Documentation #10561 available at https://www.bio-rad-antibodies.com/uploads/MSDS/10561.pdf</p>
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (H/L) (STAR207...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M284078:160321'

Printed on 26 May 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)