

Datasheet: VPA00017

Description:	RABBIT ANTI AIF
Specificity:	AIF
Format:	Purified
Product Type:	PrecisionAb™ Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Rat, Mouse N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by ion exchange chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	Peptide corresponding to amino acids 517 - 531 of human AIF
External Database Links	UniProt: O95831 Related reagents Entrez Gene: 9131 AIFM1 Related reagents
Synonyms	AIF, PDCD8

Specificity **Rabbit anti Human AIF antibody** recognizes human apoptosis-inducing factor (AIF), also known as human apoptosis-inducing factor 1, mitochondrial and programmed cell death protein 8. AIF is a 613 amino acid member of the FAD-dependent oxidoreductase family. Rabbit anti Human AIF antibody reacts with an internal region of AIF (within the region aa 517-531). AIF is widely expressed in human tissues and localizes in the mitochondria. AIF translocates to the nucleus upon induction of apoptosis. Multiple isoforms of AIF are generated by alternative splicing. Rabbit anti Human AIF antibody recognizes isoforms 1, 2, 3 and 5. It is not expected to recognize isoform 4 or 6 as they lack the immunogen sequence.

Mutations in the AIFM1 gene can result in various disease conditions including [COXPD6](#) and [Cowchock syndrome](#).

Rabbit anti Human AIF antibody recognizes human AIF as a single band of ~67 kDa in multiple human cell line lysates and cross reacts with rodent cell line lysates tested.

Western Blotting	Anti AIF antibody detects a band of approximately 67 kDa in HEK293 cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M284826:160321'

Printed on 01 May 2018