

## Datasheet: VPA00380

|                      |                         |
|----------------------|-------------------------|
| <b>Description:</b>  | RABBIT ANTI EIF3B       |
| <b>Specificity:</b>  | EIF3B                   |
| <b>Format:</b>       | Purified                |
| <b>Product Type:</b> | PrecisionAb™ Polyclonal |
| <b>Isotype:</b>      | Polyclonal IgG          |
| <b>Quantity:</b>     | 100 µl                  |

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

|                  | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| Western Blotting | ■   |    |                | 1/1000             |

**PrecisionAb antibodies have been extensively [validated for the western blot application](#).** The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

|                                 |                                                                                                                                                                      |
|---------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Target Species</b>           | Human                                                                                                                                                                |
| <b>Species Cross Reactivity</b> | Reacts with: Mouse<br><b>N.B.</b> Antibody reactivity and working conditions may vary between species.                                                               |
| <b>Product Form</b>             | Purified IgG - liquid                                                                                                                                                |
| <b>Preparation</b>              | Rabbit polyclonal antibody purified by affinity chromatography                                                                                                       |
| <b>Buffer Solution</b>          | Phosphate buffered saline                                                                                                                                            |
| <b>Preservative Stabilisers</b> | 0.09% Sodium Azide (NaN <sub>3</sub> )<br>2% Sucrose                                                                                                                 |
| <b>Immunogen</b>                | Synthetic peptide corresponding to the C-terminal region of human EIF3B                                                                                              |
| <b>External Database Links</b>  | <b>UniProt:</b><br><a href="#">P55884</a> <a href="#">Related reagents</a><br><br><b>Entrez Gene:</b><br><a href="#">8662</a> EIF3B <a href="#">Related reagents</a> |
| <b>Synonyms</b>                 | EIF3S9                                                                                                                                                               |

|                                      |                                                                                                                                                                                                                                                                                                                                                |
|--------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Specificity</b>                   | <p><b>Rabbit anti Human EIF3B antibody</b> recognizes the EIF3B, also known as eIF3 p110, eIF3 p116, eukaryotic translation initiation factor 3, subunit 9, hPrt1 and PRT1.</p> <p>Rabbit anti Human EIF3B antibody detects a band of 116 kDa. This antibody has been extensively validated for western blotting using whole cell lysates.</p> |
| <b>Western Blotting</b>              | Anti EIF3B detects a band of approximately 116 kDa in NIH/3T3 cell lysate                                                                                                                                                                                                                                                                      |
| <b>Instructions For Use</b>          | Please refer to the <a href="#">PrecisionAb western blotting protocol</a> . For additional information on secondary antibody dilution and exposure time see product web page.                                                                                                                                                                  |
| <b>Storage</b>                       | Store undiluted at -20°C, avoiding repeated freeze thaw cycles                                                                                                                                                                                                                                                                                 |
| <b>Shelf Life</b>                    | As supplied, 12 months from date of despatch                                                                                                                                                                                                                                                                                                   |
| <b>Acknowledgements</b>              | PrecisionAb™ is a trademark of Bio-Rad Laboratories                                                                                                                                                                                                                                                                                            |
| <b>Health And Safety Information</b> | Material Safety Datasheet documentation #10045 available at:<br>Antibody (10045): <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10045.pdf</a>                                                                                                                            |
| <b>Regulatory</b>                    | For research purposes only                                                                                                                                                                                                                                                                                                                     |

## Related Products

### Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

|                                  |                                                                                                                                         |                  |                                                                                                                                                 |               |                                                                                                                                                     |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>North &amp; South America</b> | Tel: +1 800 265 7376<br>Fax: +1 919 878 3751<br>Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a> | <b>Worldwide</b> | Tel: +44 (0)1865 852 700<br>Fax: +44 (0)1865 852 739<br>Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a> | <b>Europe</b> | Tel: +49 (0) 89 8090 95 21<br>Fax: +49 (0) 89 8090 95 50<br>Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a> |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|

'M286908:160718'

**Printed on 01 May 2018**