

Datasheet: VPA00085

Description:	RABBIT ANTI FLIP GAMMA/DELTA
Specificity:	FLIP GAMMA/DELTA
Other names:	I-FLICE
Format:	Purified
Product Type:	PrecisionAb™ Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	A peptide corresponding to aa191-209 at the C-terminus of human FLIP.
External Database Links	UniProt: O15519 Related reagents Entrez Gene: 8837 CFLAR Related reagents
Synonyms	CASH, CASP8AP1, CLARP, MRIT
Specificity	Rabbit anit Human FLIP gamma/delta antibody recognizes the apoptosis related protein

designated FLIP (FLICE inhibitory protein) or I-FLICE.

FLIP interacts with the adaptor protein FADD, and with caspase-8 (FLICE) and caspase-10 (FLICE2) to inhibit apoptosis induced by death receptors including CD95, DR3, TRAIL-R and TNFR1.

FLICE inhibitory proteins were initially identified in virus and human and were accordingly designated as v-FLIPs and c-FLIPs. Four splice variants of c-FLIP have been identified and are classified as FLIP alpha, beta, delta and gamma. This antibody recognizes the human delta and gamma variants but not the alpha or beta variants.

FLIP is widely expressed in human tissues. Rabbit anti Human FLIP gamma/delta antibody detects FLIP as a band of 25 kDa in human Jurkat, Ramos and Raji cell line lysates.

Western Blotting	Rabbit anti FLIP gamma/delta detects a band of approximately 25 kDa in Jurkat cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M284857:160321'

Printed on 02 May 2018