

Datasheet: VPA00031

Description:	GOAT ANTI GNB2L1
Specificity:	GNB2L1
Format:	Purified
Product Type:	PrecisionAb™ Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Goat polyclonal antibody by purified ammonium sulfate precipitation followed by antigen affinity chromatography
Buffer Solution	TRIS buffered saline
Preservative Stabilisers	0.02% Sodium Azide (NaN ₃) 0.5% Bovine Serum Albumin
Immunogen	Peptide sequence C-NTLGVCKYTVQDESH from an internal region of GNB2L1
External Database Links	UniProt: P63244 Related reagents Entrez Gene: 10399 GNB2L1 Related reagents
Specificity	Goat anti Human GNB2L1 antibody recognizes human guanine nucleotide-binding protein subunit

beta-2-like 1 (GNB2L1), also known as cell proliferation-inducing gene 21 protein, human lung cancer oncogene 7 protein, receptor for activated C kinase or RACK1. GNB2L1 is a 317 amino acid member of the WD repeat G protein beta family, containing two [WD](#) repeats. GNB2L1 acts as a substrate for Src protein-tyrosine kinase ([Chang et al. 2002](#)) and as a regulator of signaling pathways ([Sato et al. 2014](#), [Long et al. 2014](#)).

Goat anti Human GNB2L1 antibody recognizes human GNB2L1 as a band of ~35 kDa in multiple human cell lysates by western blotting under reducing conditions. The antibody also cross-reacts with rodent cell line lysates tested.

Western Blotting	Anti GNB2L1 antibody detects a band of approximately 35 kDa in HEK293 cell lysates.
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10058 available at: Antibody (10058): https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (H/L) (5160-2104...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M304766:170714'

Printed on 30 Apr 2018