

Datasheet: VPA00059

Description:	GOAT ANTI GRP94
Specificity:	GRP94
Format:	Purified
Product Type:	PrecisionAb™ Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Goat polyclonal antibody purified by affinity chromatography
Buffer Solution	TRIS buffered saline
Preservative	0.02% Sodium Azide (NaN ₃)
Stabilisers	0.5% Bovine Serum Albumin
Immunogen	Peptide with sequence C-KEGVKFDSEKTKE from the internal region of the protein sequence.
External Database Links	<p>UniProt: P14625 Related reagents</p> <p>Entrez Gene: 7184 HSP90B1 Related reagents</p>
Synonyms	GRP94, TRA1

Specificity **Goat anti Human GRP94 antibody** recognizes GRP94 a 94 kDa glucose regulated protein (GRP94), also known as tumor rejection antigen 1 (TRA1), heat shock protein 90 kDa beta member 1 (HSP90B1) and endoplasmic reticulum resident molecular chaperone, which is a paralog of cytoplasmic hsp90. GRP94 assists in the folding of membrane-bound proteins including the toll-like receptors, certain subclasses of integrins and some secreted proteins such as IgG.

Goat anti Human GRP94 detects GRP94 as a band of 100 kDa in a wide range of human cell line lysates tested. Cross reactivity with rodent cell lines is also observed.

Western Blotting Goat anti GRP94 detects a band of approximately 100 kDa in MCF7 cell lysates.

Instructions For Use Please refer to the [PrecisionAb western blotting protocol](#). For additional information on secondary antibody dilution and exposure time see product web page.

Storage Store undiluted at -20°C, avoiding repeated freeze thaw cycles.

Shelf Life As supplied, 12 months from date of despatch.

Acknowledgements PrecisionAb™ is a trademark of Bio-Rad Laboratories.

Health And Safety Information Material Safety Datasheet documentation #10058 available at: Antibody (10058): <https://www.bio-rad-antibodies.com/uploads/MSDS/10058.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Donkey Anti Sheep IgG (STAR88...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M284840:160321'

Printed on 04 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)