

Datasheet: VPA00229

Description:	RABBIT ANTI THYMOPOETIN
Specificity:	THYMOPOETIN
Format:	Purified
Product Type:	PrecisionAb™ Polyclonal
Isotype:	Polyclonal IgG
Quantity:	100 µl

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	▪			1/1000

PrecisionAb antibodies have been extensively [validated for the western blot application](#). The antibody has been validated at the suggested dilution. Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Further optimization may be required dependant on sample type.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Rabbit polyclonal antibody purified by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Immunogen	KLH conjugated synthetic peptide between 243-269 amino acids from the central region of human thymopoetin
External Database Links	<p>UniProt: P42167 Related reagents</p> <p>Entrez Gene: 7112 TMPO Related reagents</p>
Synonyms	LAP2

Specificity	<p>Rabbit anti Human thymopoietin antibody recognizes the beta isoform of thymopoietin, also known as lamina-associated polypeptide 2, isoform beta, TP beta and thymopoietin-related peptide isoform beta.</p> <p>Encoded by the TMPO gene, thymopoietins reside in the nucleus and may play a role in the assembly of the nuclear lamina, helping to maintain the structural organization of the nuclear envelope. Mutations in TMPO are associated with dilated cardiomyopathy and this gene has been found up regulated in several types of cancer. Alternatively spliced transcript variants encoding different isoforms have been noted for TMPO (provided by RefSeq, May 2010).</p> <p>Rabbit anti Human thymopoietin antibody detects a band of 51 kDa. The antibody has been extensively validated for western blotting using whole cell lysates.</p>
Western Blotting	Anti thymopoetin detects a band of approximately 51 kDa in HEK293 cell lysates
Instructions For Use	Please refer to the PrecisionAb western blotting protocol . For additional information on secondary antibody dilution and exposure time see product web page.
Storage	Store undiluted at -20°C, avoiding repeated freeze thaw cycles.
Shelf Life	As supplied, 12 months from date of despatch.
Acknowledgements	PrecisionAb™ is a trademark of Bio-Rad Laboratories.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: Antibody (10040): https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Rabbit IgG (H/L) (STAR208...) [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M284933:160321'

Printed on 01 May 2018